


# AMERITUS ENTRAL™ FEEDING TUBES

FOR ENTERAL FEEDING



Neonatal & Infant  
PRODUCTS

## ENTERAL FEEDING TUBES

Ameritus offers a full line of dedicated enteral feeding tubes in various materials, french sizes and lengths, optimizing patient safety, ease of use and versatility.

# AMERITUS ENTRAL™ FEEDING TUBES

## FOR ENTERAL FEEDING

### ASSURE PATIENT SAFETY

#### Features

The Ameritus ENTRAL™ Feeding Tubes have been designed to provide optimum patient safety while improving the standard of care and ease of use for the clinicians

- The connection on the feeding tube is specifically designed to help prevent dangerous misconnections to IV lines and Luer Syringes. The oral hub is compatible with standard oral tip syringes and enteral only extension sets (Please refer to the complete Ameritus enteral line of oral syringes and extension sets)
- Orange color coding provides easy visual recognition of the enteral connection
- Radiopaque stripe for X-ray visualization
- Ergonomic and smooth ridges on the lower section of the connector enhance grip for effective connection and disconnection of the feeding tube
- French size, length, and **lot number** clearly printed on the feeding tube and centimeter markings are printed **every centimeter**
- Unique, open-ended, rounded tip with one large eyelet for optimal flow
- Ribbed plug for connector to ensure a sealed closure
- Sterile, single packed
- All materials are Latex Free, DEHP Free, and BPA Free

#### ORDERING INFORMATION - POLYURETHANE

Catalog Number	Description	Quantity
AEG-40P-40	<b>4.0 FR, 40 cm (16"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-60P-40	<b>6.0 FR, 40 cm (16"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-80P-40	<b>8.0 FR, 40 cm (16"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-40P-50	<b>4.0 FR, 50 cm (20"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-60P-50	<b>6.0 FR, 50 cm (20"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-80P-50	<b>8.0 FR, 50 cm (20"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-40P-80	<b>4.0 FR, 80 cm (32"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-60P-80	<b>6.0 FR, 80 cm (32"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-80P-80	<b>8.0 FR, 80 cm (32"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-100P-80	<b>10.0 FR, 80 cm (32"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-120P-80	<b>12.0 FR, 80 cm (32"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-100P-110	<b>10.0 FR, 110 cm (43"), Polyurethane</b> Enteral Feeding Tube	10/bx
AEG-120P-110	<b>12.0 FR, 110 cm (43"), Polyurethane</b> Enteral Feeding Tube	10/bx

\*For Neonatal and Pediatric Use Only

#### ORDERING INFORMATION - SILICONE

Catalog Number	Description	Quantity
AEG-40S-40	<b>4.0 FR, 40 cm (16"), Silicone</b> Enteral Feeding Tube	10/bx
AEG-60S-40	<b>6.0 FR, 40 cm (16"), Silicone</b> Enteral Feeding Tube	10/bx
AEG-80S-40	<b>8.0 FR, 40 cm (16"), Silicone</b> Enteral Feeding Tube	10/bx
AEG-40S-50	<b>4.0 FR, 50 cm (20"), Silicone</b> Enteral Feeding Tube	10/bx
AEG-60S-50	<b>6.0 FR, 50 cm (20"), Silicone</b> Enteral Feeding Tube	10/bx
AEG-80S-50	<b>8.0 FR, 50 cm (20"), Silicone</b> Enteral Feeding Tube	10/bx
AEG-40S-80	<b>4.0 FR, 80 cm (32"), Silicone</b> Enteral Feeding Tube	10/bx
AEG-60S-80	<b>6.0 FR, 80 cm (32"), Silicone</b> Enteral Feeding Tube	10/bx
AEG-80S-80	<b>8.0 FR, 80 cm (32"), Silicone</b> Enteral Feeding Tube	10/bx

\*For Neonatal and Pediatric Use Only

#### ORDERING INFORMATION - PVC

Catalog Number	Description	Quantity
AEG-50V-40	<b>5.0 FR, 40 cm (16"), PVC</b> Enteral Feeding Tube	50/bx
AEG-65V-40	<b>6.5 FR, 40 cm (16"), PVC</b> Enteral Feeding Tube	50/bx
AEG-80V-40	<b>8.0 FR, 40 cm (16"), PVC</b> Enteral Feeding Tube	50/bx

Please contact our customer service department for additional information

**AMERITUS**  
Medical Products

Manufactured for and Distributed by:

Kentec Medical, Inc.


17871 Fitch, Irvine, CA 92614

TEL 949.863.0810 FAX 949.833.9730


TOLL FREE (USA only) 800.825.5996

sales@ameritusmedical.com

www.ameritusmedical.com



Neonatal & Infant  
PRODUCTS



**POLYURETHANE** Enteral Feeding Tube  
**AEG-xxP-xx**  
\*For Neonatal and Pediatric Use Only



**SILICONE** Enteral Feeding Tube  
**AEG-xxS-xx**  
\*For Neonatal and Pediatric Use Only



**PVC** Enteral Feeding Tube  
**AEG-xxV-xx**

DISTRIBUTED BY: